

Main: (800) 561-6070 Direct: (325) 356-2233



Ref. # 21257

1 of 1



League City- PowerLoad- Quote

Quote Number: 11190414

Remit to:

Stryker Sales, LLC
21343 NETWORK PLACE
CHICAGO IL 60673-1213
USA

Version: 1

Prepared For: LEAGUE CITY EMS

Attn:

Rep:

Samir Sehic

Email:

samir.sehic@stryker.com

Phone Number:

Mobile:

+1 9402977015

Quote Date: 10/01/2025

Expiration Date: 12/30/2025

Contract Start: 10/01/2025

Contract End: 09/30/2026

Delivery Address		Sold To - Shipping		Bill To Account	
Name:	LEAGUE CITY EMS	Name:	LEAGUE CITY EMS	Name:	CITY OF LEAGUE CITY
Account #:	20206214	Account #:	20206214	Account #:	20127363
Address:	144 PARK AVE STE 100	Address:	144 PARK AVE STE 100	Address:	
	LEAGUE CITY		LEAGUE CITY		
	Texas 77573-2446		Texas 77573-2446		

Equipment Products:

#	Product	Description	Qty	Sell Price	Total
1.0	639005550001	MTS POWER LOAD	1	\$35,557.00	\$35,557.00
Equipment Total:					\$35,557.00

ProCare Products:

#	Product	Description	Qty	Sell Price	Total
2.0	SYK-FS-PL-INSTALL	SYK SVC TECH POWERLOAD INSTALL 10/02/2025 - 10/01/2026	1	\$385.00	\$385.00
ProCare Total:					\$385.00

Price Totals:

Estimated Sales Tax (0.000%):	\$0.00
Shipping and Handling:	\$0.00
Grand Total:	\$35,942.00

Prices: In effect for 30 days

Terms: Net 30 Days



League City- PowerLoad- Quote

Quote Number: 11190414

Remit to:

Stryker Sales, LLC
21343 NETWORK PLACE
CHICAGO IL 60673-1213
USA

Version: 1

Prepared For: LEAGUE CITY EMS

Attn:

Rep:

Samir Sehic

Email:

samir.sehic@stryker.com

Phone Number:

Mobile:

+1 9402977015

Quote Date: 10/01/2025

Expiration Date: 12/30/2025

Contract Start: 10/01/2025

Contract End: 09/30/2026

Shipping & Handling Includes:

Standard freight, special packaging, semi rigging cranes, labor & delivery of equipment to final location, removal of all packaging, pre-delivery site check, education/training

Terms and Conditions:

Deal Consummation: This is a quote and not a commitment. This quote is subject to final credit, pricing, and documentation approval. Legal documentation must be signed before your equipment can be delivered. Documentation will be provided upon completion of our review process and your selection of a payment schedule. Confidentiality Notice: Recipient will not disclose to any third party the terms of this quote or any other information, including any pricing or discounts, offered to be provided by Stryker to Recipient in connection with this quote, without Stryker's prior written approval, except as may be requested by law or by lawful order of any applicable government agency. A copy of Stryker Medical's terms and conditions can be found at https://techweb.stryker.com/Terms_Conditions/index.html.

Zac Jordan – Marketing Manager
Stryker Medical
3800 E. Centre Ave. Portage, MI 49002



Date: December 29, 2025

Re: Power-LOAD Cot Fastener Sole Source Information

To Whom It May Concern,

Stryker Medical certifies that we are the sole manufacturer of the Stryker EMS Power-LOAD cot fastener (Model 6390). This correspondence is to inform you of the characteristics of the Power-LOAD Cot Fastener. These characteristics can be broken down into two primary categories: Independent Qualification and Ease of Use.

The Stryker EMS Power-LOAD (Model 6390) cot fastening system is mounted within the patient compartment and is intended to aid in the loading/unloading of patients. The Stryker Power-LOAD is the only powered cot fastening system that meets the following:

Independent Qualification

- IPX6: The system is tested to withstand powerful water jets.
- IEC 60601-1 and IEC 60601-1-2: The Power-LOAD cot fastener conforms to industry standards for mechanical and electrical safety for medical electrical devices, as well as electromagnetic compatibility and immunity.
- BS EN-1789:2007, clause 4.5.9 when used with Power-PRO cot and X-Restraints: This is a European dynamic crash test which subjects a 50th percentile dummy to nominal 10g deceleration for a minimum of 50ms. Following the test there shall be no sharp edges or danger to the safety of persons in the road ambulance.
- SAE J3027 compliant when used with a Stryker Power-PRO cot and X-Restraints

Ease of Use

- Device must provide a linear guide when loading and unloading the cot.
- Device must allow for remote actuation from Power-PRO foot end controls.
- Device must engage to the cot during loading and unloading, providing a means of lifting and lowering.
- Device must have a safe working load of 870 lbs. and be capable of lifting patients weighing up to 700 lbs.
- Device must be mounted inside the patient compartment to prevent environmental exposure and corrosion.
- Device must be power washable.
- Device must be capable of inductively charging the Stryker SMRT cot battery.

Please forward any further questions to your Stryker sales representative.

Sincerely,

A handwritten signature in black ink, appearing to read "Zac Jordan".

Zac Jordan
Marketing Manager



League City- Stretcher- Quote

Quote Number: 11221376

Remit to:

Stryker Sales, LLC
21343 NETWORK PLACE
CHICAGO IL 60673-1213
USA

Version: 1

Division:

Medical

Prepared For: LEAGUE CITY EMS

Rep:

Samir Sehic

Attn:

Email:

samir.sehic@stryker.com

Phone Number:

Mobile:

+1 9402977015

Quote Date: 12/01/2025

Expiration Date: 12/31/2025

Delivery Address		Sold To - Shipping		Bill To Account	
Name:	LEAGUE CITY EMS	Name:	LEAGUE CITY EMS	Name:	CITY OF LEAGUE CITY
Account #:	20206214	Account #:	20206214	Account #:	20127363
Address:	144 PARK AVE STE 100	Address:	144 PARK AVE STE 100	Address:	300 W WALKER ST
	LEAGUE CITY		LEAGUE CITY		LEAGUE CITY
	Texas 77573-2446		Texas 77573-2446		Texas 77573-3837

Equipment Products:

#	Product	Description	Qty	Sell Price	Total
1.0	650705550001	6507 POWER PRO 2, HIGH CONFIG	1	\$39,314.00	\$39,314.00
2.0	650707000002	KIT, ALVARIUM BATTERY, SERVICE	1	\$794.25	\$794.25
Equipment Total:					\$40,108.25

Price Totals:

Estimated Sales Tax (0.000%):	\$0.00
Shipping and Handling:	\$523.20
Grand Total:	\$40,631.45

Prices: In effect for 30 days

Terms: Net 30 Days



League City- Stretcher- Quote

Quote Number: 11221376

Remit to:

Stryker Sales, LLC
21343 NETWORK PLACE
CHICAGO IL 60673-1213
USA

Version: 1

Division:

Medical

Prepared For: LEAGUE CITY EMS

Rep:

Samir Sehic

Attn:

Email:

samir.sehic@stryker.com

Phone Number:

Mobile:

+1 9402977015

Quote Date: 12/01/2025

Expiration Date: 12/31/2025

Shipping & Handling Includes:

Standard freight, special packaging, semi rigging cranes, labor & delivery of equipment to final location, removal of all packaging, pre-delivery site check, education/training

Terms and Conditions:

Deal Consummation: This is a quote and not a commitment. This quote is subject to final credit, pricing, and documentation approval. Legal documentation must be signed before your equipment can be delivered. Documentation will be provided upon completion of our review process and your selection of a payment schedule. Confidentiality Notice: Recipient will not disclose to any third party the terms of this quote or any other information, including any pricing or discounts, offered to be provided by Stryker to Recipient in connection with this quote, without Stryker's prior written approval, except as may be requested by law or by lawful order of any applicable government agency. A copy of Stryker Medical's terms and conditions can be found at https://techweb.stryker.com/Terms_Conditions/index.html.

Zac Jordan – Marketing Manager
Stryker Medical
3800 E. Centre Ave. Portage, MI 49002



Date: December 29, 2025

Re: Power-PRO Ambulance Cot Sole Source Information

To Whom It May Concern,

Stryker Medical certifies that we are the sole manufacturer of the Stryker EMS Power-PRO XT ambulance cot (Model 6506). This correspondence is to inform you of the characteristics of the Power-PRO ambulance cot. These characteristics can be broken down into two primary categories: Independent Qualification and Ease of Use and Maintenance

Independent Qualification

- IPX6: The system is tested to withstand powerful water jets.
- IEC 60601-1 and IEC 60601-1-2: The Power-PRO XT cot conforms to industry standards for mechanical and electrical safety for medical electrical devices, as well as electromagnetic compatibility and immunity.
- BS EN-1789:2007, clause 4.5.9 when used with Power-LOAD cot fastener and X-Restraints: This is a European dynamic crash test which subjects a 50th percentile dummy to nominal 10g deceleration for a minimum of 50ms. Following the test there shall be no sharp edges or danger to the safety of persons in the road ambulance.
- Meets SAE J3027 dynamic crash test safety standards when used with Power-LOAD/Performance-LOAD cot fastener and X-Restraints.

Ease of Use and Maintenance

- The cot has a weight capacity of 700 lbs.
- When unloading with the manual release handle, the cot utilizes hydraulic dampening. Thus, the cot will not abruptly jar the operator or patient.
- The battery is placed at the foot-end of the stretcher.
- The cot legs power-retract in 2.4 seconds which speeds load time.
- The cot provides a load height of 36" and is operator adjustable to match the deck height of individual ambulances.
- The foot-end of the cot provides lifting bars and operator controls at two different heights, thus providing optimum ergonomics to most operator heights.
- The foot-end of the cot contains a large battery indicator light which displays amber or green which indicates battery level. A warning is given by a flashing amber light, providing the operator the time to change the battery before full depletion of power.
- The model 6506 has 6" x 2" sealed casters and bearings.
- The cot features a foot-end mounted hourly usage meter. This tool can be used to determine the timing of preventative maintenance checks.
- The cot features powder-coating of the entire aluminum frame (including the patient handling surfaces), thus eliminating aluminum oxidation throughout the cot.
- All caster bearings are sealed.
- The cot is power-washable.

Please contact your Sales Representative for further information.

Sincerely,

A handwritten signature in blue ink, appearing to read "Zac Jordan".

Marketing Manager